

# TESTIMONIALS

## WHAT PEOPLE ARE SAYING ABOUT OUR COOL FLOW OR® PRODUCTS

*"The surgeons just absolutely love the cooling vests. We've had to order one for each because they were fighting over who gets to use it. Great product, works wonderful and they love it!"*

- Terre Haute Regional Hospital in Indiana

*"I am very pleased with my Polar Surgeon's [fitted] cooling vest. It keeps me cool during surgery at my clinic. The unit is compact and works just as good or better than the huge \$2600 unit that they have available for me to use at the hospital."* - Dr. Hornstein in New Jersey

*"We love the new nylon cooling vest! We really like the new material that allows the vest to be wiped down. We will be ordering more."* - Wellstar Health in Georgia

*"The Dr. is very happy with his Cooling Vest and now we have another surgeon asking for his own vest."* - Cathy in Florida

*"I love my Cool Flow OR® Vest. It's a game changer for me. I recommend it to all my colleagues who work with lead vests because it makes the day so much easier and faster." - Dr. Raj in Washington Courthouse, Ohio*

Read more customer reviews at [polarproducts.com](http://polarproducts.com).



Polar Products is a family-owned company with over 40 years of experience manufacturing effective, high-quality and affordable hot/cold therapy and body cooling products.

[polarproducts.com](http://polarproducts.com) | 800.763.8423

## NON-TETHERED COOLING OPTIONS

### COST-EFFECTIVE COOLING FOR THE ENTIRE O.R. STAFF!

Discreet pockets hold your choice of our reusable cooling pack technologies:

#### KOOL MAX® WATER-BASED COOLING PACKS

require a freezer and offer the highest level of cooling in any environment for up to 3-4 hours.




#### COOL58® PHASE CHANGE PACKS

freeze at 58°F and cool at a constant 58°F for up to 3 hours. Ideal choice when freezers are not available; can be charged in ice water or refrigerator.




Vests are available in cotton or easy-to-clean nylon.



**ADJUSTABLE ZIPPER COOLING VEST**



**FASHION COOLING VEST**



**FASHION COOLING SCARF**

# COOL FLOW O.R.®

## COOLING VESTS & ACCESSORIES FOR SURGEONS & STAFF



## THE COOL FLOW O.R.® DIFFERENCE

### KEEP YOUR OR STAFF COOL



#### LIGHTWEIGHT, DURABLE VESTS

Provide maximum cooling coverage and comfort.

#### SAFE, WATER-TIGHT COOLING RESERVOIR

Features a neoprene gasket seal and stainless steel compression latches to ensure the floor of your operating room stays dry.

#### VERSATILE SYSTEM

Two vest styles are available. The fitted and adjustable vests both come in cotton. The adjustable vest is also available in engineered breathable nylon that is easy to clean and can be wiped down.

#### TEMPERATURE FLOW CONTROL

The digital control allows you to modify the flow of the water, modifying the temperature for extra customization and comfort.

#### SAFE OPERATION

The system runs off 12 volt DC power.

SEE OUR COMPLETE PRODUCT LINE AT [POLARPRODUCTS.COM](http://POLARPRODUCTS.COM)



# COOL FLOW O.R.

The lightweight Cool Flow O.R.® Circulating Water System keeps surgeons and staff cool and comfortable by continuously circulating cold water through tubing sewn within the vest or accessory.

## SURGICAL COOLING VEST SYSTEM FEATURES

### SAFE, WATER-TIGHT COOLING RESERVOIR:

Neoprene gasket seal, high density insulation and stainless steel compression latches ensure the system is water-tight.

### QUICK-CONNECT COUPLINGS:

Easily detach the vest from the cooling reservoir to allow the surgeon to have freedom of movement as needed.

### ATTACHED WHEELS & RETRACTABLE HANDLE:

Easy to move around for safety and convenience.

### INSULATED WATER LINES:

Detachable, insulated cool water lines from the reservoir to vest are available in any length. (10-foot length is standard with 5-foot extensions)



## ALSO AVAILABLE

### WHEELED BASKET CART :

Ideal for easy mobility, hanging the vest(s) and storing the cooling reservoir and tubing when not in use.



### ARCTIC CHILLER SYSTEM:

Eliminates the need for ice and allows you to set the desired temperature.



## COOLING RESERVOIRS

We offer three sizes of cooling reservoirs to fit your needs. The cooling reservoirs are constructed with high-performance coolers, and available in 13, 19 and 30 quart. The larger the reservoir, the more ice you can use and therefore the longer your cooling will last.



## FITTED VEST

The fitted style vest is constructed of a high-quality fabric with clean lines and a small amount of stretch which allows for a comfortable, slim fit. Sizes range from S to 4XL.

97% cotton / 3% spandex blend (only option available)



## ADJUSTABLE VEST

The adjustable style vest will fit a larger range of people and features chest and waist adjustment straps to cinch the vest to the body. Use the velcro straps at the waist to modify to the specific size. Sizes range from S to 3XL.

Available in 100% cotton or engineered breathable nylon, which is easy-to-clean and can be wiped down.



## MULTI-PERSON SYSTEMS AVAILABLE

**COOL TWO, FOUR OR SIX PEOPLE AT THE SAME TIME WITH A SINGLE COOLING RESERVOIR!**

**CONTACT US TO START YOUR FREE 45 DAY TRIAL**

[polarproducts.com](http://polarproducts.com) | 800.763.8423 | [polar@polarproducts.com](mailto:polar@polarproducts.com) | GSA Contract #36F79723D0207