



# POLAR Products

**FREE 45 Day Trial Program for CoolOR® Vest Systems**

**Thank you for your interest!**

If you have any questions, contact Jacob Graessle: jacob@polarproducts.com or 1.800.763.8423. Please email or fax (330.253.4233) the completed form.

**Trial Contact:** \_\_\_\_\_

**Hospital:** \_\_\_\_\_

**Street Address:** \_\_\_\_\_

**City, State, Zip:** \_\_\_\_\_

**Phone:** \_\_\_\_\_

**E-mail:** \_\_\_\_\_

**A hard copy of a no charge purchase order is required to initiate a free trial.**



**Complete CoolOR® Cooling Vest System includes:**

- Cooling Vest (choose style & size below)
- Water-Tight Cooling Reservoir (choose style & size below) featuring a neoprene gasket seal, high density insulation and stainless-steel compression latches
- Drain Accessory
- Digital Temperature Flow Control
- 12 Volt/3 Amp Grounded Transformer
- Lithium-ion Battery Pack with Charger
- Quick Disconnect Dry Couplings
- 8 Foot Detached Insulated Tubing Extension



**CHOOSE COOLING RESERVOIR STYLE AND SIZE:**



**19 Quart (#CORS-19)**

- \_\_\_ with integrated wheels and retractable handle
- \_\_\_ Wheeled Basket Cart instead of integrated wheels and handle



**30 Quart (#CORS-30)**

- \_\_\_ with integrated wheels and retractable handle
- \_\_\_ Wheeled Basket Cart instead of integrated wheels and handle



**13 Quart (#CORS-13)**

with carrying handle (no wheels or cart)

**CHOOSE COOLING VEST STYLE, MATERIAL AND SIZE:**



**Adjustable Vest (#CFA)**

- Size:**
- \_\_\_ S/M      \_\_\_ XL/2XL
  - \_\_\_ M/L      \_\_\_ 2XL/3XL
  - \_\_\_ L/XL

- Material:**
- \_\_\_ Breatheable, Easy-Clean Nylon
  - \_\_\_ Cotton (97%) / Spandex (3%)



**Fitted Vest (#CFF)**

- Size:**
- \_\_\_ S      \_\_\_ XL      \_\_\_ 4XL
  - \_\_\_ M      \_\_\_ 2XL
  - \_\_\_ L      \_\_\_ 3XL

- Material:**
- Cotton (97%) / Spandex (3%) (only one choice available)

**PLEASE READ AND SIGN TERMS AND CONDITIONS ON REVERSE**



## TERMS AND CONDITIONS OF FREE 45 DAY TRIAL:

To initiate the FREE 45 day trial, we will need both the no charge PO and this completed trial form (*signed below*). We require a hard copy of the no charge PO to be sent via e-mail or fax. We will send out a brand new vest and system for the trial.

At the conclusion of the trial:

1. If you decide to keep the system, we will bill you for it.
2. If you decide to return the system, it must be returned as soon as possible at your cost. There will be no charges from Polar Products.

Please e-mail or fax completed form and no charge PO to:

Jacob Graessle

E-mail: [jacob@polarproducts.com](mailto:jacob@polarproducts.com)

Fax: 330.253.4233 • Phone: 800.763.8423

I agree to provide performance feedback within 45 days of receipt of this system, unless a special request is submitted in writing to Polar Products to extend this time if needed. I agree to either authorize payment for this system at or before the end of the 45 day period or return the system to Polar at my cost. If the product is damaged or lost while in my possession for the trial period, the trialing facility will be responsible for the cost of the product. I understand that failure to comply with these conditions will result in the automatic billing for the system at 60 days on the submitted purchase order and payment will be due immediately.

*Signature of Authorized Person Requesting Free Trial:*

Printed Name: \_\_\_\_\_ Title: \_\_\_\_\_

Signature: \_\_\_\_\_ Date: \_\_\_\_\_

*Signature of Authorized Purchasing Person:*

Printed Name: \_\_\_\_\_ Title: \_\_\_\_\_

Signature: \_\_\_\_\_ Date: \_\_\_\_\_